

Eco Serie II Therapy Chair

Economic and extremely stable, the Actualway[®] Eco Série II Treatment Chair has been specifically developed to provide comfort, security and wellness during long periods of treatment such as dialysis, chemotherapy, blood transfusions and other medical treatments.

Battery

ABS trav

Crash button

Foot rest motor

Coloured castor caps

FEATURES

- Easy to use
- Foam thickness: 12cm
- Variable height from 53 78cm
- 3 motors: Backrest, elevation & legrest
- Adjustable head pillow
- Handset
- Trendelenburg position: -12º operated via handset
- Manual footrest adjustment

OPTIONS

Extra comfort upholstery Coated fabric: M2 or M1 Memory foam Movement handle Paper roll holder IV serum pole Electric height adjustment from 53 – 78cm

- Backrest adjustrmentLegrest adjustment
- Vertical and horizontal arm
- adjustment

 Epoxy metal frame
- 4 Castors with 10cm diameter

0000."

- Easy clean
- Weight capacity: 250kg

Serie II Therapy Chair

Eco Série II

Width: 55cm

Overall Width: 88cm Height: 53 – 78cm

Frame: 80 x 19 x 117cm

Overall Length: 198 – 215cm

Economic and extremely stable, the **Actualway®** Série II Treatment Chair has been specifically developed to provide comfort, security and wellness during long periods of treatment such as dialysis, chemotherapy, oncology and blood transfusions.

FEATURES

- Trendelenburg position: -12° in just 12 seconds
- Emergency button
- Adjustable height from 55 75cm
- Adjustable cervical cushion
- Easy to use handset
- 4 Castors with central brake
- Rotative castors 10cm diameter
- Manual footrest adjustment (3 motor option)

OPTIONS

Built in weight scale
Fold-away tray
IV serum pole
Movement handle
Paper roll holder
Foot rest motor
ABS base protector

For fabric options for both the Eco Série II and Série II please see the chart overleaf.

- Electrical footrest
- adjustment (4 motor option)Vertical and lateral arm rest inclination
- Stainless steel armrest support

Custom size (217 x 60cm)

Memory foam Custom foam filling Nappa coated fabric Reading lamp MCL Control box AEB Exercise bike

- Battery included
- Max load: 250kg + optional accessories

Série II

Width: 54cm Overall Width: 88cm Height: 55 – 75cm Overall Length: 217cm Frame: 65 x 25 x 105cm

option)



Electric Version Serie IV Blood Donation/ Collection Chair

The Actualway[®] Série IV Blood Donation/ Collection Chair has been specifically developed for use in a laboratories, clinics, blood services and small services. This versatile chair is also suitable for gynaecological examinations and recovery.

• Fixed front castors

components and

accessories

• 10cm Diameter castors

• High quality electrical

• Comfortable upholstery

• Class M2 fireproof and

fabric, with medical

properties efficient against bacteria (MRSA)

• Easy to clean with soap

• Wide range of colours

and water

with appropriate density

waterproof French coated

FEATURES

- Easy to use and maitain
- Quick delivery and efficient after sales service
- 3 Motors: backrest, height and legrest regulation
- Adjustable head cushion
- Easy to use handset
- Adjustable, rotating and removable arm rests (can be used as leg holders)
- Excellent stability
- Tested and approved load capacity of 200kg + accessories
- Rotating rear castors with individual locking mechanism

OPTIONS

Button to achieve horizontal position, max -4°
Stainless steel IV serum pole with 4 hooks
Stainless steel paper roll holder
Stainless steel tray with support
Memory foam
Seamless upholstery
Class M1 fireproof and waterproof fabrics
Back handle for easy movement
Special arm rest support regulation and fixation
110v Control box
Battery

Fabric options: Eco Série II, Série II and Série IV



Série IV

Min height: 50cm Max height: 70cm Overall width: 88cm Overall length: 175cm Seat rest width: 54cm Seat rest: 48 x 54 x 10cm Back rest: 76 x 54 x 10cm Leg rest: 45 x 54 x 10cm Arm rest: 48 x 15 x 4.5cm Base metal frame: 95 x 62 x 25cm

Technical specificatiuon: Eco Série II, Série II and Série IV

- Power supply: 230v/50-60Hz
- Protection classification: IP44
- Handset protection classification: IP66



For more information and for a competitive quote, please contact:

UK Distributor LINC Medical Systems LTD 7 Kingsley Street Leicgester LE2 6DY UK

Telephone +44 (0)1572 717515 **Email** sales@linc-medical.co.uk

Website www.linc-medical.co.uk

CE