

Enuresis Alarms and Sensors

LINC Code	NPC Code	Product Name	Product Photo
CMAC-001	GPB047	<p>DRI Sleeper excel bedwetting alarm</p> <ul style="list-style-type: none"> • body-worn alarm • maximum wetness detection with 9 sensing strips on both sides of sensor • responds quickly to bedwetting as the sensor • is placed inside the underwear to capture the first drops or urine • can be used with pull-ups and diapers • the plastic sensor is safe, easy to clean and dry for immediate re-use • alarm set to ring at the resonant frequency of the ear, which ear is most alerted to 	
CMAC-004	GPB073	<p>DRI Sleeper® excel Urosensor™</p> <ul style="list-style-type: none"> • replacement Urosensor 	
CMAC-003	GPB048	<p>DRI Sleeper eclipse bedwetting alarm</p> <ul style="list-style-type: none"> • wireless alarm • added comfort for the wearer as the sensor is also a transmitter- no need for a transmitter box attached to the underwear • responds quickly to bedwetting as the sensor is placed inside the underwear to capture the first drops or urine • can be used with pull-ups and diapers • the plastic sensor is safe and easy to clean and dry for immediate re-use • alarm set to ring at the resonant frequency of the ear, which the ear is most alerted to • the sensor can also be used with a second alarm for the parent's room to assist the child when they wet the bed. 	
CMAC-002	GPB072	<p>DRI Sleeper® eclipse Urosensor™</p> <ul style="list-style-type: none"> • replacement Urosensor 	

Ordering:

1. Via NHSSC: my.supplychain.nhs.uk
2. email sales@linc-medical.co.uk or tel: +44 (0) 1572 717515